

Comfortland

Pneumatic Armband

BEFORE USING THIS DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THIS DEVICE.

INDICATIONS

The Comfortland Pneumatic Armband provides stability and support for those suffering from tennis and golfer's elbow. Effective in reducing impact on the elbow joint during movement to reduce pain.

CAUTION

The application instructions listed above are intended for the use of qualified professionals only. It is not recommended for patients to attempt to apply or adjust this device unless otherwise instructed to do so by an orthopedic professional. Any questions regarding further adjustment or frequency/ duration of wear should be directed to your orthopedic professional, as they know the individual condition for which this device was prescribed.

CLEANING INSTRUCTIONS

Remove the air cell and hand wash using cold water and mild detergent. Air dry. Do not dry clean, iron, or bleach any part of this brace, as it could cause damage and void the manufacturer's warranty.

NOTICE: While every effort has been made to obtain maximum strength, function, and durability, there is no guarantee that injury will be prevented through the use of this device.

Comfortland Medical, Inc.
709 AO Smith Rd.
Mebane, NC 27302
T 866.277.3135 | F 866.277.4412
www.comfortlandmed.com

PREPARATION

1. Loosen the closure strap on the armband. Make sure that the strap stays looped through the d-ring closure for easier application (Fig. A).



APPLICATION

1. Place the armband on the arm so that the Comfortland logo is facing out. Depending on the condition and area of pain, the armband can be placed above or below the elbow (Fig. B).



2. While making sure the band stays centered, tighten the closure strap so that the band fits snug and no movement is experienced (Fig. C)

