



Product: Elastic Wrist Guard Support
Model: 2900

- Unique 3" wide elastic material adheres to hook closure strap for exact fit and easy adjustment
- Universal design fits a range of sizes
- Fits either left or right wrist



Order #
2900

Size
Universal

SIZING:

Universal design fits a range of sizes. Fits right or left.

HOW TO APPLY:

Wrap around the wrist until snug. The support should fit snug but not too tight as to be uncomfortable or restrict circulation. Attach the hook closure to secure in place.

CARE INSTRUCTIONS:

Hand wash in cold water with mild soap. Wrap in a towel to remove excess water (do not twist or wring). Air dry away from heat. Do not tumble dry. Do not bleach.

CAUTION:

Wear as recommended. Always consult a physician when encountering persistent pain, discomfort, inflammation, or swelling.

SUGG HCPCS CODE: A4466