




**Sittin Pretty Bimini Blue - 751200**

Re-Order #	Description	Size	Unit
751010	Gel-Foam Vinyl Cushion w/LSI Cover	16"W x 16"D x 2.5"H	EA
751011	Gel-Foam Vinyl Cushion w/LSI Cover	18"W x 16"D x 2.5"H	EA
751013	Gel-Foam Vinyl Cushion w/LSI Cover	18"W x 18"D x 2.5"H	EA
751012	Gel-Foam Vinyl Cushion w/LSI Cover	20"W x 16"D x 2.5"H	EA
751014	Gel-Foam Vinyl Cushion w/LSI Cover	20"W x 18"D x 2.5"H	EA
751018	Gel-Foam Vinyl Cushion Only	16"W x 16"D x 2.5"H	EA
751019	Gel-Foam Vinyl Cushion Only	18"W x 16"D x 2.5"H	EA
751021	Gel-Foam Vinyl Cushion	18"W x 18"D x 2.5"H	EA
751020	Gel-Foam Vinyl Cushion Only	20"W x 16"D x 2.5"H	EA
751039	Gel-Foam Thin-Line Vinyl Cushion w/LSI Cover	16"W x 16"D x 1.5"H	EA
751040	Gel-Foam Thin-Line Vinyl Cushion w/LSI Cover	18"W x 16"D x 1.5"H	EA
751044	Gel-Foam Thin-Line Vinyl Cushion w/LSI Cover	18"W x 18"D x 1.5"H	EA
751049	Gel-Foam Thin-Line Vinyl Cushion w/LSI Cover	20"W x 18"D x 1.5"H	EA
751042	Gel-Foam Thin-Line Vinyl Cushion w/LSI Cover	20"W x 16"D x 1.5"H	EA
914481	Gel-Foam Ultra-Thin Vinyl Cushion w/LSI Cover	16"W x 16"D x 1"H	EA
914482	Gel-Foam Ultra-Thin Vinyl Cushion w/LSI Cover	18"W x 16"D x 1"H	EA
914483	Gel-Foam Ultra-Thin Vinyl Cushion w/LSI Cover	20"W x 16"D x 1"H	EA
751025	Gel-Foam Vinyl Cushion w/Coccyx Cutout & LSI Cover	18"W x 16"D x 2.5"H	EA
754528	Gel-Foam Vinyl Cushion w/Coccyx Cutout & LSI Cover	18"W x 16"D x 3.5"H	EA
751038	Gel-Foam Thin-Line Vinyl Cushion Only	16"W x 16"D x 1.5"H	EA
751041	Gel-Foam Thin-Line Vinyl Cushion Only	18"W x 16"D x 1.5"H	EA
751045	Gel-Foam Thin-Line Vinyl Cushion Only	18"W x 18"D x 1.5"H	EA
751043	Gel-Foam Thin-Line Vinyl Cushion Only	20"W x 16"D x 1.5"H	EA
751200	Sittin Pretty Bimini Blue Gel-Foam Vinyl Cushion	16"W x 16"D x 2.5"H	EA
751201	Sittin Pretty Bimini Blue Gel-Foam Vinyl Cushion	18"W x 16"D x 2.5"H	EA
751202	Sittin Pretty Bimini Blue Gel-Foam Vinyl Cushion	18"W x 18"D x 2.5"H	EA



### Purpose

Skil-Care's Gel-Foam Cushion is a combination gel and foam cushion designed to relieve pressure. The water based gel surface provides cooling action that reduces skin-damaging perspiration.

### Important Safety Reminder

The Low-Shear cover has a safety strap that prevents the cushion from slipping off the seat. This strap must be secured to the wheelchair frame.

### Easy Care

The vinyl cover may be wiped clean using a damp cloth and spray-type cleanser. A mild disinfectant may be used. The Low-Shear outer cover can be laundered and dried at temperatures that do not exceed 180°F (82°C). Do not use bleach.