


GEL OVERLAY

Reorder #	Description	Size
914429	X-Gel Overlay	16"x16"x1"
914430	X-Gel Overlay	18"x16" x1"
781066	Low Shear II Cover for 914429, 16" x 16"	16"x16"
781065	Low Shear II Cover for 914430, 18" x 16"	18"x16"

Purpose

This pad is filled with a water-based gel and designed to be placed over foam a cushion to provide pressure-relieving flotation. The gel's cooling action reduces skin-damaging perspiration.



Important Safety Reminder

The Gel Overlay has safety ties that must be secured to the wheelchair to keep it from slipping off the foam cushion.

Easy care

The Gel Overlay wipes clean with a damp cloth and spray-type cleanser. It may be disinfected with a non-chlorine disinfectant.

Note

This package contains only the Gel Overlay. It does not include a foam cushion.