

TELEMETRY UNIT POUCHES

Reorder #	Material	Previous #	Unit
<input type="checkbox"/> 901130	Cloth	0102	36/Cs.
<input type="checkbox"/> 901131	Foam	0103	36/Cs.

These comfortable-to-wear pouches protect costly telemetry units and assist in providing more accurate readings.

EASY TO USE

Secure the top tie behind the patient's neck, and the bottom tie around the lower part of the patient's chest as shown in the illustration. For convenience in tying and for patient comfort, cut off extra tape. Place telemetry unit in pouch and close by pressing hook-and-loop tabs together.

