


# Chair Protection Pad




Reorder #	Description	Length	Unit
<input type="checkbox"/> 909260	Cushion Pad Protector	18" x 16"	12 PK
<input type="checkbox"/> 909270	Cushion Pad Protector	18" x 16"	Each

## Purpose

Prevents urine damage to lounge chairs, bedside chairs, and geri-chairs. Three-layer design soaks up and holds urine and prevents it from seeping through to the chair. Side and rear barriers provide added protection.

## Easy to Use

Secure back barrier to side barriers using the Velcro® tabs. (See illustration at right.) Place on chair. Note tag reading: ABSORBENT SIDE-PLACE FACE UP.



## Laundering

Skil-Care Chair Protection Pads may be laundered and dried using standard white goods procedures. Do not wash or dry at temperatures above 180°F (82°C). The use of bleach is not recommended.