



# Endeavor Universal Wrist/ Thumb Splint

HCPC L3807/ L3809 Approved



## Indications:

Carpal tunnel syndrome

Rheumatoid arthritis

Osteoarthritis

De Quervain's Disease

Fracture of the radius

Smith's fracture

Fracture of the shaft

Sprains and strains of the  
wrist and hand

Contact Us:

Toll Free: (866) 277-3135

Fax: (866) 277-4412

[www.comfortlandmed.com](http://www.comfortlandmed.com)



### Lightweight Design

Weighing under 4 oz., the wrist/ thumb splint provides incredible support while remaining light-weight and low-profile, helping patients get back to normal activity quickly.



### Easy to Use/ Adjust

The wrist/ thumb splint features a hook strap and D-ring closure system, allowing for easy application and adjustment without requiring tools or further assistance.



### Breathable Construction

Made of breathable soft goods, the wrist/ thumb splint promotes air-flow to the wrist and hand, speeding up the healing process as well as to help keep patients cool and dry, even after long periods of use.