

REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET TI/HI

## PATIENT CARE

**Sani-Hands® Instant Hand Sanitizing Wipes**

Contains 70% Ethyl Alcohol, Aloe &amp; Vitamin E emollients, and is clinically proven to kill 99.99% of germs.

	REORDER NO.	WIPE SIZE	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI
Medium Canister	P13472	6 x 7.5"	12/135s	18.45 lbs	0.99 ft	8/5
Large Canister	P15984	6 x 7.5"	6/220s	15.25 lbs	0.80 ft	10/5
Individual Packet	D43600	5 x 8"	10/100s	9.86 lbs	0.65 ft	9/7
Bedside Pack	P71520	5.5 x 8.4"	48/20s	11.20 lbs	0.46 ft	15/6

**Hygea® Multi-Purpose Washcloths**

Strong, disposable washcloths that are ideal for cleansing and refreshing patients.

	REORDER NO.	WIPE SIZE	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI
Premium - Solo® Softpak	J14143	11.5 x 9.5"	6/60s	11.00 lbs	0.67 ft	9/7
Solo Softpak	J14108	12 x 8"	9/60s	13.30 lbs	0.67 ft	10/5
Solo Softpak	J21396	12 x 8"	6/96s	13.30 lbs	0.72 ft	10/6
Resealable Softpak	J22750	12 x 8"	12/48s	13.30 lbs	0.67 ft	10/6
Bulk Individual Packet	U12095	11.5 x 10"	1/400s	16.70 lbs	1.56 ft	6/4

**Hygea® Flushable Personal Cleansing Cloths**The ideal flushable pre-moistened washcloth for personal general skin cleansing and incontinence care. Safe for use with sewer and well-maintained septic systems.<sup>a,b,c</sup>

	REORDER NO.	WIPE SIZE	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI
Solo Softpak	A500F48	5.3 x 6.8"	12/48s	6.89 lbs	0.37 ft	15/7

**Hygea® BZK Antiseptic Towelette (alcohol free)**

For use wherever antiseptic cleansing/refreshing is needed. Contains 0.13% Benzalkonium Chloride.

	REORDER NO.	WIPE SIZE	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI
Individual Packet	D35185	7 x 5"	20/100s	17.30 lbs	1.26 ft	9/3

**Hygea® Obstetrical Towelette**

For personal hygiene at every napkin change and whenever cleansing/refreshing is needed. Saturated with Benzalkonium Chloride, Alcohol 5%, Chlorothymol.

	REORDER NO.	WIPE SIZE	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI
Individual Packet	D74800	7.75 x 5"	10/100s	9.67 lbs	0.65 ft	9/7

<sup>a</sup>When used as directed<sup>b</sup>Based on the testing guidelines to determine flushability set by the Association of the Nonwoven Fabrics Industry (INDA)<sup>c</sup>Not for use with basement pump systems. Only flush one wipe at a time