



## **Caig Laboratories Manufacturer Warranty**

### **WARRANTY:**

CAIG will warranty its products for a period of 1 year from invoice date. Storage conditions must be adhered to.

\*\* Storage Life: If storage directions are adhered to, products will last considerably longer. Values are conservative. If coloration or separation of products occur, we recommend discarding and purchasing new product.

Storage Instructions: Store in a cool, dry place. Keep container tightly closed. Do not subject to excessive heat or cold.

Avoid Contaminating Solutions: When using bulk liquids we recommend transferring small amounts into a smaller container for use. This will avoid contamination of the larger container from dirty cleaning swabs, cloths, brushes, etc.

### **CONTACT:**

Parts Express

Ph: 1-800-338-0531 ext. 780

### **RETURN PROCESS:**

Call Parts Express to obtain RA number. Pack merchandise carefully to avoid further damage, and complete Return Form on the bottom of the invoice. Record customer and RA number on the outside of the box. Ship merchandise back to Parts Express prepaid.